



iM3s
Vital Signs Monitor

Pocket spot checks at all times

iM3s

Vital Signs Monitor



About Edan

Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
- Ultrasound Imaging
- *In-Vitro* Diagnostics
- Patient Monitoring
- Point-of-Care Testing
- Veterinary
- OB/GYN

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.



Global Headquarters :
Edan Instruments, Inc. | No.15 Jinhui Rd., Jinsha Community,
Kengzi Subdistrict, Pingshan District, Shenzhen | 518122 P.R. China
+86.755.26898326 | www.edan.com | info@edan.com.cn

U.S. and Canada inquiries :
EDAN Diagnostics, Inc. | 9918 Via Pasar, San Diego, CA 92126
+1.858.750.3066 | www.edandiagnostics.com | edan-info@edandiagnostics.com

© Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission. Not all products or features are available in all countries, contact Edan for local availability.

ENG-PM-iM3s-V1.0-20191108
For Medical use only




Small outside, smart inside

Being the smallest vital signs monitor among its kind, the iM3s presents a unique pocket-size design, delivering extraordinary portability throughout spot checks. Though small in size, it integrates all necessary vital measurements together, including SpO₂, NIBP, and optional temperature.

A Dedicated Spot Monitor

With its pleasing size, the iM3s works and operates just like a smart phone, bringing comfortable and familiar user experience to daily spot check applications. For each test, all measured parameters are displayed together on the same screen, while the test results being saved automatically in its built-in non-volatile memory.



300g



5 inch screen




IP44



Shock resistance



Wi-Fi




20,000 sets



AA Batteries

5-6 hours



Tabletop Expansion

Plugged onto an optional tabletop dock, the iM3s easily fits into regular tabletop applications, enabling more possibilities, including AC charging, optional thermal recording, optional quick temperature measurements, etc.



Connect Your Medical Records

The Wi-Fi capability allows the iM3s to connect to the MFM-CMS for centralized data management, or to communicate with HIS/EMR systems on a bi-directional basis.

